


「Index」



第9章 目录	9-1
→ 数字	9-2
→ 英文字母	9-2
→ 中文	9-6



INDEX

INDEX

数字

10mM dNTPs Mixture	2-39
100mM dNTPs (dATP,dCTP,dGTP,dTTP)	2-39
2mM dNTPs Mixture	2-39
2mM dNTPs Mixture (dATP,dCTP,dGTP,dUTP)	2-39
50×ROX reference dye	1-19

A

anti-Taq high	2-37
---------------	------

B

Blend Taq -Plus- (带 dNTPs)	2-29
Blend Taq (带 dNTPs)	2-29
Buffer	
KOD FX 用 2×PCR Buffer	2-39
KOD FX Neo 用 2×PCR Buffer	2-39
KOD -Plus- 用 buffer	2-39
KOD -Plus- 用 25mM Mg ₂ SO ₄	2-39
KOD -Plus- Neo 用 buffer	2-39
KOD -Plus- Ver.2 用 buffer	2-39

C

Can Get Signal Immunoreaction Enhancer Solution	6-6
Can Get Signal Immunoreaction Enhancer Solution 1	6-6
Can Get Signal Immunoreaction Enhancer Solution 2	6-6
Can Get Signal Immunostain Immunoreaction Enhancer Solution A	6-8
Can Get Signal Immunostain Immunoreaction Enhancer Solution B	6-8
Can Get Signal Immunostain Immunoreaction Enhancer Solution Start Set	6-8
cDNA 合成	
ReverTra Ace	3-14
ReverTra Ace qPCR RT Kit	3-10
ReverTra Ace qPCR RT Master Mix	3-8
ReverTra Ace qPCR RT Master Mix with gDNA remover	3-8
ReverTra Ace -α-	3-12
SuperPrep® II Cell Lysis & RT Kit for qPCR	3-5
SuperPrep® II Cell Lysis Kit for qPCR	3-5

D

Deoxyribonucleotides	
10mM dNTPs Mixture	2-39
100mM dNTPs (dATP,dCTP,dGTP,dTTP)	2-39
100mM dUTP	2-39
2mM dNTPs Mixture	2-39
2mM dNTPs Mixture (dATP,dCTP,dGTP,dUTP)	2-39
DNA Polymerase	
Blend Taq	2-29
Blend Taq -Plus-	2-29
KOD Dash	2-30
KOD FX	2-24
KOD FX Neo	2-22
KOD One PCR Master Mix	2-14
KOD One PCR Master Mix -Blue-	2-14
KOD -Multi & Epi-	2-26
KOD -Plus-	2-20
KOD -Plus- Neo	2-18
KOD -Plus-Ver.2	2-20
Quick Taq HS DyeMix	2-31
rTaq DNA Polymerase	2-33
rTth DNA Polymerase	2-34
DNA 纯化试剂盒	
MagExtractor -Genome-	5-4
MagExtractor -PCR & Gel Clean up-	5-10
MagExtractor -Plant Genome-	5-5
MagExtractor -Plasmid-	5-9
DNA 聚合酶	
Blend Taq	2-29
Blend Taq -Plus-	2-29
KOD Dash	2-30
KOD FX	2-24
KOD FX Neo	2-22
KOD One PCR Master Mix	2-14
KOD One PCR Master Mix -Blue-	2-14
KOD -Multi & Epi-	2-26
KOD -Plus-	2-20
KOD -Plus- Neo	2-18

KOD -Plus- Ver.2	2-20	KOD -Multi & Epi-	2-26
Quick Taq HS DyeMix	2-31	KOD -Plus-	2-20
rTaq DNA Polymerase	2-33	KOD -Plus- Mutagenesis Kit	4-8
rTth DNA Polymerase	2-34	KOD -Plus- Neo	2-18
		KOD -Plus- Neo 用 Buffer	2-39
dNTPs		KOD -Plus- Ver.2	2-20
10mM dNTPs Mixture	2-39	KOD -Plus- Ver.2 用 Buffer	2-39
100mM dNTPs (dATP,dCTP,dGTP,dTTP)	2-39	KOD -Plus- 用 25mM Mg ₂ SO ₄	2-39
100mM dUTP	2-39	KOD -Plus- 用 Buffer	2-39
2mM dNTPs Mixture	2-39	KOD SYBR® qPCR Mix	1-13
2mM dNTPs Mixture (dATP,dCTP,dGTP,dUTP)	2-39		
		L	
G		Labeling Kit	
Gene Sequencing		ScriptMAX Thermo T7 Transcription Kit	4-13
GenNext NGS Library Prep Kit	1-24		
GenNext NGS Library Quantification Kit	1-22	Ligation Kit	
GenNext RamDA-seq Single Cell Kit	1-26	Ligation high	4-10
KOD -Multi & Epi-	2-26	Ligation high Ver.2	4-11
KOD SYBR® qPCR Mix	1-13		
NSR Primer Set for human	1-26	M	
NSR Primer Set for mouse	1-26	MagExtractor -Genome-	5-4
RamDA Cell Lysis Kit	1-26	MagExtractor -PCR & Gel Clean up-	5-10
		MagExtractor -Plant Genome-	5-5
Genome 纯化系列		MagExtractor -Plasmid-	5-9
MagExtractor -Genome-	5-4	MagExtractor -RNA-	5-6
MagExtractor -Plant Genome-	5-5	MagExtractor -Viral RNA-	5-8
		Magical Trapper	5-11
H			
Handy Pestle	5-11	Magnetic Beads Kit	
		MagExtractor -Genome-	5-4
Hot Start TTx (DNA) Kit	2-35	MagExtractor -PCR & Gel Clean up-	5-10
Hot Start TTx (RNA) Kit	2-35	MagExtractor -Plant Genome-	5-5
		MagExtractor -Plasmid-	5-9
K		MagExtractor -RNA-	5-6
KOD Dash	2-30	MagExtractor -Viral RNA-	5-8
KOD FX	2-24	Magical Trapper	5-10
KOD FX Neo	2-22		
KOD FX Neo 用 2×PCR Buffer	2-39	Master Mix	
KOD FX 用 2×PCR Buffer	2-39	KOD One PCR Master Mix	2-14
KOD One PCR Master Mix	2-14	KOD One PCR Master Mix -Blue-	2-14
KOD One PCR Master Mix -Blue-	2-14	KOD SYBR® qPCR Mix	1-13

Realtime PCR Master Mix	1-21
RNA-direct Realtime PCR Master Mix	1-12
RT-PCR Quick Master Mix	2-36
THUNDERBIRD qPCR Mix	1-19
THUNDERBIRD Next qPCR Mix	1-17
THUNDERBIRD Probe One -step qRT-PCR Kit	1-8
Quick Taq HS DyeMix	2-31

mRNA 表达分析	
Hot Start TTx (DNA) Kit	2-35
Hot Start TTx (RNA) Kit	2-35
KOD SYBR® qPCR Mix	1-13
ReverTra Ace qPCR RT Kit	3-10
ReverTra Ace qPCR RT Master Mix	3-8
ReverTra Ace qPCR RT Master Mix with gDNA remover	3-8
Realtime PCR Master Mix	1-21
RNA-direct Realtime PCR Master Mix	1-12
RT-RamDA cDNA Synthesis Kit	1-26
SuperPrep® II Cell Lysis & RT Kit for qPCR	3-5
SuperPrep® II Cell Lysis Kit for qPCR	3-5
THUNDERBIRD qPCR Mix	1-19
THUNDERBIRD Next qPCR Mix	1-17
THUNDERBIRD Probe One-step qRT-PCR Kit	1-8

KOD -Plus- Mutagenesis Kit	4-8
ScriptMAX Thermo T7 Transcription Kit	4-13

Mg ²⁺	
KOD -Plus- 用 25mM Mg ₂ SO ₄	2-39

N	
NGS	
GenNext NGS Library Prep Kit	1-24
GenNext NGS Library Quantification Kit	1-22
GenNext RamDA-seq Single Cell Kit	1-26
KOD -Multi & Epi-	2-26
KOD SYBR® qPCR Mix	1-13
NSR Primer Set for human	1-26
NSR Primer Set for mouse	1-26
RamDA Cell Lysis Kit	1-26

NTP	
10mM dNTPs Mixture	2-39
100mM dNTPs (dATP,dCTP,dGTP,dTTP)	2-39
100mM dUTP	2-39
2mM dNTPs Mixture	2-39
2mM dNTPs Mixture (dATP,dCTP,dGTP,dUTP)	2-39

Nucleotides	
10mM dNTPs Mixture	2-39
100mM dNTPs (dATP,dCTP,dGTP,dTTP)	2-39
100mM dUTP	2-39
2mM dNTPs Mixture	2-39
2mM dNTPs Mixture (dATP,dCTP,dGTP,dUTP)	2-39

O	
oligo(dT) ₂₀	3-12

P	
PCR 定量	
Hot Start TTx (DNA) Kit	2-35
Hot Start TTx (RNA) Kit	2-35
KOD SYBR® qPCR Mix	1-13
ReverTra Ace qPCR RT Kit	3-10
ReverTra Ace qPCR RT Master Mix	3-8
ReverTra Ace qPCR RT Master Mix with gDNA remover	3-8
Realtime PCR Master Mix	1-21
RNA-direct Realtime PCR Master Mix	1-12
RT-RamDA cDNA Synthesis Kit	1-26
SuperPrep® II Cell Lysis & RT Kit for qPCR	3-5
SuperPrep® II Cell Lysis Kit for qPCR	3-5
THUNDERBIRD qPCR Mix	1-19
THUNDERBIRD Next qPCR Mix	1-17
THUNDERBIRD Probe One-step qRT-PCR Kit	1-8

PCR 酶	
Blend Taq	2-29
Blend Taq -Plus-	2-29
KOD Dash	2-30
KOD FX	2-24
KOD FX Neo	2-22
KOD One PCR Master Mix	2-14

KOD One PCR Master Mix -Blue-	2-14
KOD -Multi & Epi-	2-26
KOD -Plus-	2-20
KOD -Plus- Neo	2-18
KOD -Plus-Ver.2	2-20
Quick Taq HS DyeMix	2-31
rTaq DNA Polymerase	2-33
rTth DNA Polymerase	2-34
PCR 相关试剂	
anti-Taq high	2-37
Deoxyribonucleotides	2-39
PCR 用缓冲液	2-39
RNase Inhibitor	3-16
PCR 试剂盒	
RT-PCR Quick Master Mix	2-36
PCR 用缓冲液	
KOD FX 用 2×PCR Buffer	2-39
KOD FX Neo 用 2×PCR Buffer	2-39
KOD -Plus- 用 buffer	2-39
KOD -Plus- 用 25mM Mg ₂ SO ₄	2-39
KOD -Plus- Neo 用 buffer	2-39
KOD -Plus- Ver.2 用 buffer	2-39
Polymerase	
Blend Taq	2-29
Blend Taq -Plus-	2-29
KOD Dash	2-30
KOD FX	2-24
KOD FX Neo	2-22
KOD One PCR Master Mix	2-14
KOD One PCR Master Mix -Blue-	2-14
KOD -Multi & Epi-	2-26
KOD -Plus-	2-20
KOD -Plus- Neo	2-18
KOD -Plus- Ver.2	2-20
Quick Taq HS DyeMix	2-31
rTaq DNA Polymerase	2-33
rTth DNA Polymerase	2-34

Primer	
oligo(dT) ₂₀	3-12
Random Primer(9mer)	3-12
Protein kit	
ScriptMAX Thermo T7 Transcription Kit	4-13
Q	
Quick Taq HS DyeMix	2-31
qPCR	
Hot Start TTx (DNA) Kit	2-35
Hot Start TTx (RNA) Kit	2-35
KOD SYBR® qPCR Mix	1-13
ReverTra Ace qPCR RT Kit	3-10
ReverTra Ace qPCR RT Master Mix	3-8
ReverTra Ace qPCR RT Master Mix with gDNA remover	3-8
Realtime PCR Master Mix	1-21
RNA-direct Realtime PCR Master Mix	1-12
RT-RamDA cDNA Synthesis Kit	1-26
SuperPrep® II Cell Lysis & RT Kit for qPCR	3-5
SuperPrep® II Cell Lysis Kit for qPCR	3-5
THUNDERBIRD qPCR Mix	1-19
THUNDERBIRD Next qPCR Mix	1-17
THUNDERBIRD Probe One-step qRT-PCR Kit	1-8
R	
RamDA Cell Lysis Kit	1-26
Random Primer(9mer)	3-12
Realtime PCR	
Hot Start TTx (DNA) Kit	2-35
Hot Start TTx (RNA) Kit	2-35
KOD SYBR® qPCR Mix	1-13
ReverTra Ace qPCR RT Kit	3-10
ReverTra Ace qPCR RT Master Mix	3-8
ReverTra Ace qPCR RT Master Mix with gDNA remover	3-8
Realtime PCR Master Mix	1-21
RNA-direct Realtime PCR Master Mix	1-12
RT-RamDA cDNA Synthesis Kit	1-26
SuperPrep® II Cell Lysis & RT Kit for qPCR	3-5
SuperPrep® II Cell Lysis Kit for qPCR	3-5

抽提器械		高效率 PCR 酶	
Magical Trapper	5-11	KOD FX	2-24
Handy Pestle	5-12	KOD FX Neo	2-22
抽提试剂盒		KOD One PCR Master Mix	2-14
MagExtractor -Genome-	5-4	KOD One PCR Master Mix -Blue-	2-14
MagExtractor -PCR & Gel Clean up-	5-10	KOD -Multi & Epi-	2-26
MagExtractor -Plant Genome-	5-5	核糖核酸多聚酶	
MagExtractor -Plasmid-	5-9	Thermo T7 RNA Polymerase <<TT7>>	4-15
MagExtractor -RNA-	5-6	核糖核酸酶抑制剂	
MagExtractor -Viral RNA-	5-8	RNase Inhibitor	3-16
蛋白质合成		缓冲液	
ScriptMAX Thermo T7 Transcription Kit	4-13	KOD FX 用 2×PCR Buffer	2-39
定量 PCR		KOD FX Neo 用 2×PCR Buffer	2-39
Hot Start TTx (DNA) Kit	2-35	KOD -Plus- 用 buffer	2-39
Hot Start TTx (RNA) Kit	2-35	KOD -Plus- 用 25mM Mg ₂ SO ₄	2-39
KOD SYBR® qPCR Mix	1-13	KOD -Plus- Neo 用 buffer	2-39
ReverTra Ace qPCR RT Kit	3-10	KOD -Plus- Ver.2 用 buffer	2-39
ReverTra Ace qPCR RT Master Mix	3-8	聚合酶	
ReverTra Ace qPCR RT Master Mix with gDNA remover	3-8	Blend Taq	2-29
Realtime PCR Master Mix	1-21	Blend Taq -Plus-	2-29
RNA-direct Realtime PCR Master Mix	1-12	KOD Dash	2-30
RT-RamDA cDNA Synthesis Kit	1-26	KOD FX	2-24
SuperPrep® II Cell Lysis & RT Kit for qPCR	3-5	KOD FX Neo	2-22
SuperPrep® II Cell Lysis Kit for qPCR	3-5	KOD One PCR Master Mix	2-14
THUNDERBIRD qPCR Mix	1-19	KOD One PCR Master Mix -Blue-	2-14
THUNDERBIRD Next qPCR Mix	1-17	KOD -Multi & Epi-	2-26
THUNDERBIRD Probe One-step qRT-PCR Kit	1-8	KOD -Plus-	2-20
高保真性 PCR 酶		KOD -Plus- Neo	2-18
KOD FX	2-24	KOD -Plus- Ver.2	2-20
KOD FX Neo	2-22	Quick Taq HS DyeMix	2-31
KOD One PCR Master Mix	2-14	rTaq DNA Polymerase	2-33
KOD One PCR Master Mix -Blue-	2-14	rTth DNA Polymerase	2-34
KOD -Multi & Epi-	2-26	克隆	
KOD -Plus-	2-20	TARget Clone	4-6
KOD -Plus- Neo	2-18	TARget Clone -Plus-	4-6
KOD -Plus- Ver.2	2-20	KOD -Plus- Mutagenesis Kit	4-8

Ligation high 4-10

Ligation high Ver.2 4-11

连接酶试剂盒

Ligation high 4-10

Ligation high Ver.2 4-11

免疫反应增强

Can Get Signal Immunoreaction Enhancer Solution 6-6

Can Get Signal Immunoreaction Enhancer Solution 1 6-6

Can Get Signal Immunoreaction Enhancer Solution 2 6-6

Can Get Signal Immunostain Immunoreaction Enhancer Solution A 6-8

Can Get Signal Immunostain Immunoreaction Enhancer Solution B 6-8

Can Get Signal Immunostain Immunoreaction Enhancer Solution Start Set 6-8

逆转录酶

ReverTra Ace 3-14

突变导入

KOD -Plus- Mutagenesis Kit 4-8

脱氧核糖核苷酸

dNTPs Mixture 2-39

转录试剂盒

ScriptMAX Thermo T7 Transcription Kit 4-13


大幅度提高了THUNDERBIRD®性能

高效率Realtime PCR Master Mix


THUNDERBIRD® Next SYBR® qPCR Mix

同等价格可得到效率更高的Realtime PCR!

- ▶ 扩增效率提高
- ▶ 低拷贝区的检测灵敏度提高
- ▶ 特异性进一步增强
- ▶ 可对应延伸时间为10秒的高速循环
- ▶ 制备的反应液的稳定性提高
- ▶ 防止Carryover污染引起的假阳性
(与另外销售的UNG组合使用时)
- ▶ 含有校正染料, 无需添加ROX



溶液有颜色, 所以降低了
分装误差



使用『RamDA-seq™方法』从单细胞或微量RNA开始的cDNA制备试剂盒

NGS分析用的cDNA制备试剂盒

GenNext® RamDA-seq™ Single Cell Kit

Realtime PCR分析用cDNA制备试剂盒

RT-RamDA™ cDNA Synthesis Kit

✓ 单细胞及微量RNA为模板


可使用1~100个细胞或 10pg~1ng total RNA 作为input

✓ 可制备覆盖全长RNA的cDNA

现有的只使用Oligo-dT引物的方法难以分析10kb以上的RNA，使用本产品可以制备出测序reads覆盖全长RNA的cDNA

RamDA-seq™

A公司试剂盒



✓ 除了可检测poly(A) RNA外，还可检测non-poly(A) RNA

除了Oligo-dT primer外，还使用了NSR (Not so random primer)，可检测到用现有方法无法检测的non-poly(A) RNA


✓ 相比原来技术，检测到的基因数大大提高了

可检测组蛋白RNA、lncRNA、pre-mRNA、circRNA等各种类型的RNA

使用从Mouse ES细胞提取的RNA 10pg，分别用两种试剂盒 (N=2) 制备cDNA后，进行NGS分析

cDNA Synthesis Method	RamDA-seq™		A公司试剂盒	
Number of transcripts TPM>0.1	18,984	18,465	13,500	13,721

高效PCR快速选择指南



KOD -Plus- 系列	KOD FX 系列	KOD -Multi & Epi-	Blend Taq 系列
与Taq DNA聚合酶保真性的比较 (增加倍数)			
80倍	11倍		3倍
-克隆用高保真性PCR -高速PCR (30秒/kb)	-粗样品 (血液、鼠尾、植物等) -含细胞壁的微生物 (酵母、真菌等) -难扩增的目的片段 (长度>20kb, 高GC含量)	-多重PCR (~10kb) -亚硫酸氢盐PCR (~1.5kb) -可使用含I或U的引物	-普通PCR -长PCR (>5kb) -克隆PCR
			Quick Taq HS DyeMix 1倍 -普通PCR -克隆PCR -预混型 
P 2-18	P 2-22	P 2-26	P 2-31
推荐TA克隆试剂盒: Target Clone系列 (P 3-6)			
Target Clone -Plus-	Target Clone -Plus-	Target Clone -Plus-	Target Clone